




# Medtrum EasyView Website

## User Guide

*For Patients*





Medtrum Technologies Inc.

Building 3 and Building 8, No. 200,  
Niudun Road

Shanghai 201203, China

Tel: +86-21-50274781

Fax: +86-21-50274779

[www.medtrum.com](http://www.medtrum.com)

Version: 1.5

Publication date:

11 November 2024

UG889000WW-004

348838

<b>1 Introduction .....</b>	<b>1</b>
<b>1.1 Related Documents.....</b>	<b>1</b>
<b>1.2 User Safety .....</b>	<b>1</b>
1.2.1 Intended Use.....	1
1.2.2 Intended Users.....	1
1.2.3 Contraindications.....	1
1.2.4 Performance Characteristics.....	1
1.2.5 General Warnings .....	1
<b>2 Getting Started .....</b>	<b>3</b>
<b>2.1 Accessing the Website .....</b>	<b>3</b>
<b>2.2 Sign-Up.....</b>	<b>3</b>
<b>2.3 Sign-In .....</b>	<b>6</b>
<b>3 Settings .....</b>	<b>8</b>
<b>3.1 Report Settings .....</b>	<b>8</b>
<b>3.2 Data to View .....</b>	<b>9</b>
<b>4 Monitor.....</b>	<b>11</b>
<b>4.1 CGM .....</b>	<b>11</b>
<b>4.2 Insulin Pump.....</b>	<b>14</b>
<b>5 Reports .....</b>	<b>17</b>
<b>5.1 Dashboard .....</b>	<b>17</b>
<b>5.2 Daily Summary.....</b>	<b>19</b>
<b>5.3 Day by Day Overview .....</b>	<b>21</b>
<b>5.4 Sensor Overlay.....</b>	<b>23</b>
<b>5.5 Trend Analysis .....</b>	<b>24</b>
<b>5.6 Events .....</b>	<b>27</b>

# Content

5.7 Logbook .....	29
5.8 Alert .....	31
5.9 Device Settings Review .....	31
6 Records .....	35
7 Quick Print .....	36
8 Connections, Profile and Password.....	37
8.1 Connections .....	37
8.2 Profile .....	39
8.3 Security Settings.....	39
9 Icon Table .....	41

Medtrum EasyView Website enables you to view both the real-time and uploaded data of your insulin pump and continuous glucose monitoring (CGM) system. You can also print the data report for further use.

## 1.1 Related Documents

- ✓ *Medtrum EasyTouch App User Guide*
- ✓ *Medtrum EasyPatch App User Guide*
- ✓ *TouchCare® System User Guide*
- ✓ *CGM System User Guide*

## 1.2 User Safety

### 1.2.1 Intended Use

EasyView is a data management software designed to assist patients, caregivers and healthcare professionals in the review, analysis, and evaluation of data from Medtrum devices to support diabetes management.

### 1.2.2 Intended Users

EasyView is intended for use by healthcare professionals, individuals with diabetes and their caregivers.

### 1.2.3 Contraindications

None known.

### 1.2.4 Performance Characteristics

The EasyView website maintains a minimum service uptime of 99%.

### 1.2.5 General Warnings

Warning: EasyView neither measures nor interprets the data it displays and does not provide medical advice. It should be used together with guidance from a healthcare professional specializing in diabetes management.

## Introduction

Treatment decisions should always be made in consultation with a qualified healthcare professional.

**Warning:** In the event of any malfunctions with the EasyView website, please contact Medtrum customer support or your local Medtrum distributor.

**Warning:** Report any serious incident related to the device to Medtrum customer support or your local Medtrum distributor, and the competent authority of the Member State in which you are established.

## 2.1 Accessing the Website

Before accessing EasyView Website, ensure that your operating system and web browser are both supported.

### Supported Operating Systems and Web Browsers

Operating System	Web Browser
Windows 7, Windows 8, Windows 10, and Windows 11	Google Chrome 61.0 or later, Firefox 56.0.2 or later, or Microsoft Edge 41 or later
Mac OS 10.10.5 or later	Safari 8.0.8 or later


To start EasyView Website, please visit the following web address:

<https://easyview.medtrum.eu>


## 2.2 Sign-Up

The following page shows what the EasyView sign-in page looks like. If you do not have a Medtrum account, click **Sign up now** to register an account with your e-mail address first.

# Getting Started



In the displayed sign-up dialog box, choose **Patient**, select your country/region, and enter the required information to create an account.



### Description for Sign-up Parameters

Parameter	Description
E-mail address	A valid e-mail address that will be used as your account. If you have an account of Medtrum EasyTouch, EasyPatch, or EasySense app, you can directly use this account to sign in to EasyView Website.
Verification code	A 6-digit verification code sent to your e-mail address after you click <b>Get code</b> .
Full name	Your name. You are advised to enter your full name so that your healthcare providers can easily recognize you.
Password	Password of your account. Your password should contain upper case letters, lower case letters, numbers (0-9), and special characters or non-alphanumeric characters such as !, ", £, \$, %, &, *, and @.


**Note:** If there is no verification code sent to your e-mail address, please check your junk mails.

Agree to the terms of use and privacy policy of Medtrum after you have read them, and then click **Submit** to register your account.

# Getting Started

## 2.3 Sign-In

If you have registered your account in the previous step or you have an account registered on your Medtrum EasyTouch, EasyPatch, or EasySense app, you can directly sign in to EasyView Website with it.



1. Click the flag in the upper right corner of the sign-in page to choose the country/region you used upon registration.


# Getting Started



2. Enter your e-mail address used upon registration and password.
3. Click **Sign in**.

# Settings

If the sign-in is successful, you can find the following navigation bar on the left side of the displayed page. By choosing **Settings**, you can configure the settings of reports and data to view.



## 3.1 Report Settings

Choose **Settings** > **Report Settings** to set reports.

You can toggle **Show alerts** and **Show Glucose Value** on or off to decide whether to display them in the report.

Set **Glucose Target Range**. Its value will be used in the following graphs: **Dashboard**, **Daily Summary**, **Day by Day Overview**, **Sensor Overlay**, **Trend Analysis**, and **Logbook**.

**Report Settings (mmol/L)**

Show alerts		<input type="checkbox"/> Off		Show Glucose Value		<input type="checkbox"/> Off	
Very High		13.9		High		10.0	
Low		3.9		Very Low		3.0	
Hypoglycaemia		3.0					
		Glucose Target Range		Time Period		Post Meal Analysis	
		Low      High		hours		hours	
Breakfast	Before	3.9	7.2	07:00 - 08:00		01:00 - 03:00	
	After	5.6	8.9				
Lunch	Before	3.9	7.2	12:00 - 13:00		01:00 - 03:00	
	After	5.6	8.9				
Dinner	Before	3.9	7.2	18:00 - 19:00		01:00 - 03:00	
	After	5.6	8.9				
Night-time		5.6 - 8.9		23:00 - 06:00			

**Save changes**    Restore defaults


Set **Time Period** for breakfast, lunch, dinner and nighttime. They will be used in graphs of **Day by Day Overview** and **Logbook**.

Click **Save changes** to save your settings.


## 3.2 Data to View

Choose **Settings > Data To View**.


**Data To View**



**Uploaded Data**  
The data is from EasyTouch App



**Real-time Data**  
The data is from EasySense App or EasyPatch App.

 Delete Data

If you are using a TouchCare® System with a Medtrum EasyTouch APP, please select **Uploaded Data**.


If you are using an EasySense or EasyPatch APP, please select **Real-time Data**.

**Note:** Please refer to *Medtrum EasyTouch App User Guide* to learn about how to use EasyTouch with your TouchCare® System.

## Settings


**Note:** You can click **Delete Data** to delete both real-time and uploaded data.

Choose **Monitor** from the navigation bar to display the **Monitor** page.



The Monitor page will display all your real-time device information if you are using a Medtrum product connected via Bluetooth to a Medtrum app that has Internet access.

## 4.1 CGM



If you are using a glucose sensor, you will be able to see the real-time sensor data in the **CGM** area.

# Monitor


## CGM Data Description

Data	Description
Time of the last sensor reading	Date and time of the last sensor reading sent to cloud.
Transmitter SN	The serial number of the current transmitter.
Status	<ul style="list-style-type: none"><li>● <b>Not Calibrated:</b> The sensor is not calibrated after warm-up.</li><li>● <b>Monitoring:</b> The sensor and transmitter function normally, and the transmitter is sending the sensor glucose to the app.</li><li>● <b>Calibration Needed:</b> A new calibration is required.</li><li>● <b>Calibration Error:</b> The last calibration failed.</li><li>● <b>Lost Sensor:</b> The sensor signal is lost.</li><li>● <b>Sensor Expired:</b> The sensor has been used over 14 days.</li><li>● <b>No Readings:</b> The sensor functions abnormally.</li><li>● <b>Sensor Failure:</b> The sensor has failed.</li><li>● <b>Transmitter Error:</b> The transmitter functions abnormally.</li><li>● <b>Connecting Sensor:</b> The transmitter is not properly connected with an active sensor.</li></ul>

Data	Description
	<ul style="list-style-type: none"> <li>● <b>Charging Battery:</b> The transmitter battery is being charged.</li> <li>● <b>Battery Charged:</b> The transmitter battery is fully charged.</li> </ul>
Sensor life left	The rest functioning time of the sensor.
Last sensor reading	The latest sensor reading your app received.
Transmitter battery	Power consumption of the transmitter.
Sensor glucose curve	The green curve in the graph indicating the trend of your sensor glucose.
Calibration record	The blood drip icon in the graph indicating the calibrated glucose reading.
New sensor	The green square icon in the graph indicating the use of a new sensor.
Warm up	CGM warm-up period marked in grey.
Alert	Indicates a CGM alert is triggered.

# Monitor

## 4.2 Insulin Pump



If you are using an insulin pump, you will see the real-time pump data in the **Insulin Pump** area.

### Insulin Pump Data Description


Data	Description
Time of the last pump data	Date and time of the last pump data sent to cloud.
Pump Base SN	The current serial number of the pump base.
Status	<ul style="list-style-type: none"><li>● <b>Delivering Basal:</b> The patch pump is delivering basal insulin.</li><li>● <b>Delivering Normal Bolus:</b> The patch pump is delivering a normal bolus.</li><li>● <b>Delivering Extended Bolus:</b> The patch pump is delivering an extended bolus.</li></ul>

Data	Description
	<ul style="list-style-type: none"> <li>● <b>Suspend:</b> Insulin delivery is suspended.</li> <li>● <b>Lost Pump:</b> The PDM has lost pump signal.</li> <li>● <b>Occlusion Detected:</b> The patch pump is occluded.</li> <li>● <b>Empty Reservoir:</b> No insulin left in the patch pump.</li> <li>● <b>Patch Expired:</b> The reservoir patch has expired.</li> <li>● <b>Patch Error:</b> The reservoir patch is not working properly.</li> <li>● <b>Patch Battery Depleted:</b> The reservoir patch battery is depleted.</li> <li>● <b>Pump Base Error:</b> The pump base is not working properly.</li> <li>● <b>Auto Mode:</b> The auto mode is enabled.</li> </ul>
Reservoir patch life left	The reservoir path life left.
Insulin left	The amount of insulin left in the reservoir.
Basal	The current basal rate.

## Monitor

Data	Description
Total daily basal	The total daily basal insulin that has been delivered.
Last bolus	The last bolus dose.
Active insulin	The amount of active insulin.
Total daily bolus	The total daily bolus insulin that has been delivered.


Choose **Reports** from the navigation bar to display the **Reports** page.



You can find the analysis of the sensor or pump data in the following reports: **Dashboard**, **Daily Summary**, **Day by Day Overview**, **Sensor Overlay**, **Trend Analysis**, **Events**, **Logbook**, **Alert** and **Device Settings Review**.

If you specify the ending date on one of the reports, the ending date of the other reports will be synchronized automatically.

## 5.1 Dashboard



The dashboard report is used for reviewing the statistics of CGM readings and insulin pump delivery. It provides detailed data for analyzing and managing diabetes.

The first part of the report shows the overall glucose statistics: average BG, CGM reading, daily insulin delivery, carbs-intake, and exercise time as well as standard deviation (SD) and number of input (#) of each. An estimated A1C and GMI is calculated from SG values.


# Reports

BG		CGM		Insulin		Carbs		Exercise	
Average		Average		Average Daily Dose		Average Carbs/Day		Average time/Day	
8.1 mmol/L		6.4 mmol/L		13.60 U		189.2 g		--Hours -- min	
SD	#	SD	#	SD	# days	SD	#	SD	#
2.16	57	1.83	4770	3.05	6.0	13.5	42	--	--

Estimated A1C: 5.7% calculated from SG values  
GMI (%): 6.1

Note: Estimated A1C does not replace lab measurement

The second part displays the percentile report, in which you can view the SG distribution at certain daily time points during the selected days.




The third part of the report includes detailed glucose and insulin delivery statistics. You can also view the glucose and insulin distribution of the selected days.

SG	BG/Calibration (per day)	4.7 / 2
	AUC high > 10mmol/L 10 mmol/L	0.1
	AUC low < 3.9mmol/L 3.9 mmol/L	0
	Highest value	20.4 mmol/L
	Lowest value	2.8 mmol/L
	<3.1 mmol/L	0.2%
	<3.9 mmol/L	0.8%
	3.9 mmol/L-10 mmol/L	95.9%
	>10 mmol/L	3.3%
	>13.9 mmol/L	1.2%
CV		0.31
Sensor	Sensor Usage	2/9.9days/70.7%
Pump	Daily Basal	11.56U/77.2 %
	Daily Bolus	3.41U/22.8 %
	Average Bolus No. (per day)	0
	Average Suspend time (per day)	4 min
	Average Reservoir usage (per day)	2.0 days


Very High 1.2%  
High 3.3%  
Target 95.9%  
Low 0.8%  
Very Low 0.2%

Legend: >13.9, >10, 3.9-10, <3.9, <3.1 Unstimmmol/L

Overall: Time In Range



The fourth part of the report includes Auto Mode statistics. You can view the percentage of time during which your BG stays in the target range when the Auto Mode is enabled. You can also view the amount of auto basal and auto correction bolus delivered from the table.



You can choose to view data from the past 7 days, 14 days, 30 days, 90 days, or custom duration.

Duration:

7 days
14 days
30 days
90 days
Customize

Date:

03-12-2022
-
09-12-2022

You can print this report by simply clicking **Print** in the upper right corner of the page.


## 5.2 Daily Summary

Dashboard
Daily Summary
Day by Day Overview
Sensor Overlay
Trend Analysis
Events
Logbook
Alert

Device Settings Review


# Reports

This report presents the sensor glucose curve, calibration marks, basal rates, bolus doses, alerts, carbohydrate records, target SG in auto mode, time changes made in device settings, manual insulin injections, and exercise records of a selected day.



Statistics						
	Glucose		Insulin		Events	
	BG	SG				
Avg. (mmol/L)±SD	7.2±1.33	6.6±1.19	Total daily dose (U)	16.95	Carbs (g)	40
Min (mmol/L)	5.6	4.3	Bolus (U)/%	1.10/6%	# of meals	2
Max (mmol/L)	8.5	9.3	Basal (U)/%	15.85/94%	Avg±SD	20.0±0.0
Target (%)	100	100	Food Bolus (U)	0	Total exercise time	30
Low (%)	0	0	Correction Bolus (U)	0	# of exercises	1
High (%)	0	0	Food+corr (U)	0.65	Avg. Intensity	30
# of readings	4	418	Normal Bolus (U)	0.65		
# of hypo/duration(min)	0	0/0	Extended Bolus (U)	0		
# of alerts	--	2	Temp Basal (U)	0		
			Duration of suspension(min)	0		
			# of pump alarms	0		
Auto Mode						
Time in Auto Mode(%)	23%					
Avg.(mmol/L)±SD	7.1 ± 1.16 mmol/L					
Target Range(%)	100%					
Above Range(%)	0%					
Below Range(%)	0%					
# of Auto Mode Exits	0					

By moving the cursor along the graph, you will see detailed information of each bolus delivery.



The default date is always the last date before the day you open the website.


You can change the date by clicking

Date: 06-12-2022 - 12-12-2022

and choose the date on the calendar.

You can print this report by simply clicking **Print**.

## 5.3 Day by Day Overview




This report presents the meter blood glucose readings, sensor glucose curves, insulin delivery curves, and carbs intake of 7 or 14 days. The default period is 7 days.

# Reports



You can move your mouse cursor to the blood drop icon of the BG you wish to view.




You may select to view data from the past 7 days, 14 days, or custom duration.

Duration:
☒ 7 days
☐ 14 days


Date:

You can print this report by simply clicking **Print**.

## 5.4 Sensor Overlay




This report displays the sensor data overlay and meter blood glucose readings within a given number of days (7 days at most) before a selected date. All the daily SG curves and meter blood glucose readings are displayed in an overlap graph so that you can easily see the pattern of glucose levels in a given period. If you print it, you will also see the daily average SG, maximum SG, minimum SG, and the number of hypoglycemia episodes.



By moving the cursor along the curves, you can see more detailed information.

# Reports



Although the Sensor Overlay page displays only the data of the past 7 days, you can specify the start date or the end date to view the corresponding graph.

Date: 03-12-2022 - 09-12-2022

Week: Mon Tue Wed Thu Fri Sat Sun

You can print this report by simply clicking **Print**.

## 5.5 Trend Analysis

Dashboard Daily Summary Day by Day Overview Sensor Overlay **Trend Analysis**

Events Logbook Alert Device Settings Review

This report shows the trend graph and statistics within a given number of days (90 days at most).

The Glucose Trend graph and statistics include:

### ➤ Basic statistics

- **Avg.:** The average of all the sensor readings.

- **Max:** The highest sensor reading.
- **Min:** The lowest sensor reading.
- **Med:** The Midpoint of all the sensor readings.
- **Q1:** 25% of sensor readings are lower than this value.
- **Q3:** 75% of sensor readings are lower than this value.
- **IQR:** The interquartile range (IQR) is the difference between **Q1** and **Q3**.
- **SD:** The Standard deviation (SD) measures the amount of variability or dispersion, from the individual sensor readings to the average.
- **CV:** The Coefficient of variation (CV) is a measure of relatively variability.
- **SE Mean:** The Standard Error (SE) Mean measures how spread out the sensor readings are.
- **Distribution %:** The distribution of the sensor readings. The high and low limits are set in the report settings.
- **Stability**
  - **AUC above limit:** The Area under a Curve (AUC) above limit is a measure of the high glucose above the high limit.
  - **AUC below limit:** The AUC below limit is a measure of the low glucose below the low limit.
  - **# of high / low:** Shows how many times the sensor readings are higher than the high limit or lower than the low limit.
- **BG statistics**
  - **Avg:** The average of all the BG readings.
  - **SD:** The Standard deviation (SD) measures the amount of variability or dispersion, from the individual BG reading to the average.
  - **# of Hypo:** Shows how many times the BG readings are lower than 3.1 mmol/L (56 mg/dL).
  - **BG readings:** The number of the BG readings.

# Reports

The Insulin trend graph and statistics include:


- **Avg. Daily Total:** The average of total insulin delivered per day.
- **Avg. Daily Basal:** The average of daily Basal delivered per day.
- **Avg. Daily Bolus:** The average of daily Bolus delivered per day.
- **Avg. # of bolus/Day:** The average number of times to deliver a bolus per day.
- **Avg. bolus(U)/Each:** The average insulin delivered for each Bolus.

The Carbs statistics include:

- **Carbs (g):** The carbohydrates intake
- **# of meals:** The number of meals.


The Exercise statistics include:

- **Duration (min):** The total duration of exercise
- **# of Exercises:** The number of exercises.



Statistics											
Statistics	Insulin						Carbs		Exercise		
	Avg. Daily Total (U)	Daily Basal	Daily Bolus	Avg. # of bolus/Day	Avg. bolus(U)/Each		Carbs (g)	# of meals	Duration (min)	# of Exercises	
3-12	11.85	10.5	88.6%	1.35	11.4%	0	--	0	0	0	0
4-12	0	0	--	0	--	0	--	0	0	0	0
5-12	26.75	17.35	64.9%	9.4	35.1%	5	1.9	92	5	0	0
6-12	43.15	38.25	88.6%	4.9	11.4%	1	4.9	15	1	0	0
7-12	62.9	55.4	88.1%	7.5	11.9%	0	--	0	0	0	0
8-12	54.7	52.25	95.5%	2.45	4.5%	0	--	0	0	0	0
9-12	43	42.65	99.2%	0.35	0.8%	0	--	0	0	0	0
Avg/Total	40.4	36.05	89.3%	4.35	10.7%	1	4.35	54	3	0	0

By moving the cursor along the graphs, you can see more detailed information.



You can choose to view data from the past 7 days, 14 days, 30 days, 90 days or custom duration.

Duration:

Interval:

Date:
 -











Week:
☒ Mon
☒ Tue
☒ Wed
☒ Thu
☒ Fri
☒ Sat
☒ Sun

You can print this report by simply clicking **Print**.


## 5.6 Events


# Reports


You can view events recorded on your EasyTouch, EasySense, or EasyPatch app and those recorded in your PDM if you have uploaded PDM data to the cloud.


<div>Delete</div> <div>Create</div>						
	Time	Event Type	Details	Notes	Source	Action
<input type="checkbox"/>	09-12-2022 14:48:20	BG	6.0 mmol/L		app	 
<input checked="" type="checkbox"/>	08-12-2022 10:56:54	BG	10.0 mmol/L		app	 
<input type="checkbox"/>	08-12-2022 09:43:47	BG	13.9 mmol/L		app	 
<input type="checkbox"/>	06-12-2022 15:19:06	BG	6.9 mmol/L		app	 
<input checked="" type="checkbox"/>	06-12-2022 14:04:54	BG	6.5 mmol/L		app	 

To create a new event from the website, click **Create**, record the event in the displayed dialog box, and then click **Create**.



To edit an existing event, click  in the **Action** column, edit the event in the displayed dialogue box, and then click **Save**.




To delete events from the page, select them and then click **Delete** or  in the **Action** column.

You can choose to view data from the past 7 days, 14 days, 30 days, 90 days or custom duration.


Duration:	<b>7 days</b>	14 days	30 days	90 days	Customize				
Event Filter:	<b>All</b>	Carbs	Exercise	HbA1c	Urine Ketone	BG	Injection	Health	Other
Date	23-10-2024 - 29-10-2024								

## 5.7 Logbook




This report presents the meter blood glucose values, bolus doses, pump status and other events for each hour of a specified period.

# Reports



You can move the horizontal scroll bar to see all the data throughout the 24 hours. BG readings within the normal range are marked in green, while high BG readings and low BG readings are marked in yellow and red respectively. The color for data visualization is entirely dependent on your preference of the overall glucose target. Numbers inside the rectangular box (**Manual Bolus or Modified Bolus**) represents either manual bolus or bolus given through Bolus Calculator but differs from the suggested amount.

By clicking a bolus value in the report, you will see the detailed information of this bolus delivery.



You can choose to view data from the past 7 days, 14 days or custom duration.

Duration: **7 days** 14 days

Date: 03-12-2022 - 09-12-2022

Week: ☒ Mon ☒ Tue ☒ Wed ☒ Thu ☒ Fri ☒ Sat ☒ Sun

You can print this report by simply clicking **Print**.

## 5.8 Alert

Dashboard
Daily Summary
Day by Day Overview
Sensor Overlay
Trend Analysis
Events
Logbook
Alert
Device Settings Review

This report presents all the Pump, CGM and PDM alarms and alerts within a specified period. You can choose the duration of the alert history, the alert type, and the specific alerts that you need.

Duration:
7 days
14 days
30 days
90 days
Customize

Type:
Pump
CGM
PDM

Date
28-06-2023 - 04-07-2023

Alert
Patch Expired

Time	Type	Alert
04-07-2023 02:35:12	Pump	Patch Expired
03-07-2023 18:35:21	Pump	Patch Expired
03-07-2023 17:35:21	Pump	Patch Expire In 1 Hour
03-07-2023 04:35:25	Pump	Patch Expiration Advisory

## 5.9 Device Settings Review

Dashboard
Daily Summary
Day by Day Overview
Sensor Overlay
Trend Analysis
Events
Logbook
Alert
Device Settings Review

This report presents the pump and CGM settings of a specified day. Additionally, the PDM or app settings are displayed depending on the data source selected.

# Reports



Standard													24h Total: 28.30U
Start	00:00	02:00	03:00	05:00	06:00	07:00	09:00	11:00	12:00	15:00	17:00	19:30	
U/H	1.00	1.05	1.10	1.20	1.30	1.45	1.30	1.20	0.95	1.20	1.25	1.30	
Start	21:30	23:00											
U/H	1.20	1.10											
Pattern A													24h Total: 28.15U
Start	00:00	02:00	03:00	06:30	14:00	16:30	19:00	23:00					
U/H	1.00	1.10	1.15	1.20	1.25	1.15	1.25	1.05					

**Preset temp basal**

HeavyEx	Sick	Temp1
74.00%	200.00%	4.00U/H
240	240	150

Only the basal patterns and preset temp basal that have been preprogrammed will be displayed on the **Device Settings Review** page.

# Reports

Bolus						
Maximum Bolus: 25U		Bolus Calculator: On			Active Insulin Time: 3H	
IC Ratios(g/U)						
Start	00:00	07:30	11:00	12:30	16:30	19:30
Carbs	7.8	6.7	7.2	8	7.5	7
Insulin Sensitivity mmol/L						
Start		00:00				
BG		2.8				
BG Target						
From		00:00				
Range		5.0-6.1				
Bolus Reminder (Off)						
From		--	--	--	--	--
To		--	--	--	--	--

Pump								
Pump Base SN: 3F2D4GT0			Patch Expiration: Off			Low Reservoir: Insulin 12U		
Out of Range: Off			Command Reminder: Off			Auto Off: On 12H		
Daily Max: 80U			Hourly Max: 35 U					
Alarms & Alerts on APP			Alarms on App			Alerts on App		
Audio	Vibrate	Pump Light	Repeat	Snooze	Tone	Repeat	Snooze	Tone
On	Off	On	5	60	Default	10	60	Default

CGM			
Transmitter SN: A1B2C3D4		Sensor Expiration: Off	Glucose Alerts: Off
Factory Calibration: On			
Glucose Limits			
From		00:00	
Low		4.2	
High		10.0	

# Reports

Predictive Alerts (hh:mm)			Rate Alerts (mmol/L/min)					
High	On	30	Rise	On	0.220			
Low	On	30	Fall	On	0.220			
Alerts for Low			Alerts for High			Alerts for Others		
Repeat	2		Repeat	2		Repeat	5	
Snooze	15		Snooze	30		Snooze	30	
Tone	Default		Tone	Notification		Tone	Notification	
Alerts for Hypo								
Tone	Notification							

APP		
Name: Bob Papgeorgiou	Audio Options: Audio	Glucose Unit: mmol/L
Weight Unit: kg	Height Unit: cm	
EasyLoop		
Low Suspend: Off	Predictive Low Suspend: Off	
Auto Mode		
Target SG: 5.6 mmol/L	Auto Meal Handling: On	Auto Mode: Off
Total Daily Dose: 62 U		Weight: 65 kg


If there's more than one device setting profile saved, you may select from date and time to view another profile.

Time:	12-12-2022 09:45:21
	12-12-2022 09:45:21
	09-12-2022 09:52:24

You can print this report by simply clicking **Print**.

Choose **Records** from the navigation bar to open the **Records** page.

On the **Records** page, you can find the medical records, if any, that your healthcare provider has entered.



You are not allowed to edit your medical records.

Records

Print

Profile:

Username: <a href="#">abc@gmail.com</a>	Name: <a href="#">abcd</a>	Type of diabetes: Unspecified
Height:	Weight:	Age:
Duration:		

Complications:

Concomitant diseases	Date of diagnosis	Suggestions
No Data		

Prescription:


Medicine	Dose	Notes
No Data		

Other suggestions: 09-12-2022 - 11-12-2022

HCP	Date	Suggestions
No Data		

# Quick Print

Choose **Print** from the navigation bar to display the **Quick Print** page.



On this page you can quickly print different reports. To do that, choose the duration, and the end date, specify the period, select the reports you want to print, and click **Create Report** to print.

Select duration: **One Week** One Month One Quarter

End date: **Yesterday** Last Saturday Last Sunday Customize

Period: 05-12-2022 - 11-12-2022

Week: ☒ Mon ☒ Tue ☒ Wed ☒ Thu ☒ Fri ☒ Sat ☒ Sun

Quick Print Settings

<input checked="" type="checkbox"/>	All		
<input checked="" type="checkbox"/>	Dashboard		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Daily Summary		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Day by Day Overview		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Sensor Overlay		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Trend Analysis		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Logbook		<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Event Summary	Date Item	<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Alert	PDM Pump CGM	<input checked="" type="checkbox"/> 05-12-2022 - 11-12-2022
<input checked="" type="checkbox"/>	Device Settings Review		


Create Report

# Connections, Profile and Password

Click the down arrow next to your name to display the drop-down menu.



The red dot indicates you have connection requests.




## 8.1 Connections

Choose **Connections** from the drop-down menu to display the **Connections** page.


In this page, you will find all the connection requests made by your healthcare providers or family members.

You can either agree to or deny the request.


# Connections, Profile and Password



If you choose to agree to the request, your device data will be shared with others.



If you want to delete a connection, click the small cross in the upper right corner of the connection card, and choose **Delete** on the displayed screen. After that, the person of the deleted connection card is not allowed to access to your device data.




**Note:** As a patient, your account cannot send out connection requests.

### 8.2 Profile

Choose **Profile** from the drop-down menu to display the **Profile** page. You can edit or delete your profile picture on this page.

Click the **Change** button next to your username to update it. Updating username does not cause data loss.

If your username is updated, the username displayed on your Medtrum apps will remain the old one. It is refreshed upon your next sign-in with the new account.



The screenshot displays the Medtrum Profile page. At the top, the username is shown as 'bing...' followed by '@medtrum.com' and a blue 'Change' button. To the right is a profile picture placeholder with 'Edit profile picture' and 'Delete profile' buttons. The form is divided into two columns. The left column contains: 'Full name:' with a text field containing 'MedtrumTest'; 'Country/Region:' with a dropdown menu showing 'United Kingdom (current)'; 'Height:' with a text field containing '0' and a unit dropdown set to 'cm'; 'Profession:' with a 'Select' dropdown; 'Date of birth:' with a text field and a calendar icon; 'Mobile:' with a text field; 'Email:' with a text field; and 'Address:' with a large text area. The right column contains: 'Gender:' with a 'Select' dropdown; 'Weight:' with a text field containing '0' and a unit dropdown set to 'kg'; 'Organization:' with a text field; 'Bio:' with a text field; 'Telephone:' with a text field; and 'Postcode:' with a text field. At the bottom right are 'Save changes' and 'Cancel' buttons.

### 8.3 Security Settings

Choose **Account and security** from the drop-down menu to configure security settings.

If you toggle **Multi-factor authentication (MFA)** on, a dynamic password will be required at each time when your account is performing the sign-in operation.

## Connections, Profile and Password

The current email for receiving Multi-Factor Authentication (MFA) codes is displayed, and you can click **Change** to update it if needed. Note that changing the MFA email will not affect the account name.

You can change your password in the **Change Password** area. Once you have set a new password, click **Save changes** to save the configuration.

If you want to deregister your account, click **Delete** to delete your account permanently. Once you delete your account, you will no longer be able to use this account to sign in to EasyFollow or EasyView Website.

Account and security

Multi-factor authentication (MFA)

Every time you log in, in addition to entering a password, you also need to enter a dynamic password to further enhance the security of your account.

MFA Email:  [Change](#)

☐ Off

Change Password

Old password

New password

Re-enter password

[Forgot your password?](#)









[Save changes](#) [Cancel](#)

Delete The Account Permanently

Please note that once you delete your account, all the data in your account will be lost, and you will no longer be able to use this account to log in to EasyTouch, EasySense, EasyPatch or EasyView website.

[Delete](#)

## Icon Table

	Authorized representative in the European Community
	Importer
	Unique device identifier
	Indicates the item is a Medical Device
	CE mark
	Date of manufacture
	Manufacturer
	Consult instructions for use or consult electronic instructions for use

